

Welcome!

Thank you for choosing our AED cabinet. This manual is designed to provide detailed guidance on installation, use, and maintenance of the product, ensuring a safe and reliable experience. Please read carefully to make the most of your cabinet.

Contents

1. Introduction
 2. Included Components
 3. Wall Mounting Instructions
 4. Hardware and Assembly Guide
 5. Lighting System Instructions
 6. Troubleshooting
 7. Product Specifications
 8. Care and Maintenance
 9. Frequently Asked Questions
-

1. Introduction

This AED cabinet is specifically designed to house an Automated External Defibrillator (AED) securely while providing rapid access during emergencies. Built with durable materials and practical features, the cabinet ensures long-term reliability, safety, and ease of use.

Key benefits include:

- Durable steel construction for strength and longevity.
 - Emergency lighting for visibility in critical situations.
 - Easy wall mounting with included hardware.
-

2. Included Components

Upon receiving the product, verify that the package contains all necessary components:

- AED Cabinet
 - AED Sign (for clear visibility)
 - Security Seals
 - Inspection Cards (to log maintenance checks)
 - Handle Screws (separately packed)
 - Wall Mounting Hardware:
 - Screws
 - Wall Anchors
 - Washers
 - External Lighting System (pre-installed inside)
 - AED Instruction
 - Instruction Sheet (this document)
-

3. Wall Mounting Instructions

1. **Choose an Optimal Location:** Select a wall that is sturdy, visible, and easily accessible in emergencies.
 2. **Mark Drill Points:**
 - Position the cabinet where you plan to install it and mark the mounting holes with a pencil.
 - Ensure the cabinet is level using a spirit level (a tool that indicates horizontal alignment). Alternatively, you can use a smartphone app designed for leveling.
 3. **Drill Holes:**
 - Drill holes at the marked points to fit the wall anchors provided.
 4. **Insert Wall Anchors:**
 - Push the wall anchors into the drilled holes until flush with the wall.
 5. **Attach the Cabinet:**
 - Align the cabinet with the prepared holes.
 - Use the screws and washers to secure the cabinet firmly to the wall.
 6. **Verify Stability:**
 - Test the cabinet to ensure it is securely mounted and does not move or wobble.
-

4. Hardware and Assembly Guide

- **Labeling:**
 - Each hardware bag is labeled for easy identification.
 - Refer to the diagram below to match hardware components with their corresponding parts.

[Insert diagram of labeled hardware components here]

- **Assembly Tips:**
 - Use appropriate screws for specific assembly tasks, such as attaching the handle and mounting hardware.
 - Tighten all screws and nuts securely to avoid rattling or loose parts.
-

5. Lighting System Instructions

The lighting system is pre-installed inside the cabinet to prevent damage during transportation. Once you receive and unpack the cabinet you need to unscrew the lighting system from its secure position inside the cabinet and carefully place it on the top of the cabinet.

- **Features:**
 - The cabinet is equipped with a red emergency light that automatically activates when the door is opened, ensuring visibility in low-light conditions.
- **Battery Installation:**
 - Open the cabinet and locate the battery compartment near the light system.
 - Insert 9V battery (not included) following the polarity indicators.
- **Switch Operation:**
 - The switch located on the wiring box controls the light. Toggle the switch to the "ON" position for activation.
- **Troubleshooting:**
 - If the light does not function:

- Verify batteries are correctly installed.
 - Inspect connections for loose or disconnected wires.
 - Replace the batteries if they appear drained.
-

6. Troubleshooting

- **Loose Parts:**
 - Use a standard screwdriver or wrench to tighten any loose screws or nuts.
 - **Lighting Issues:**
 - Check battery orientation and replace it as necessary.
 - Ensure the switch is toggled to the "ON" position.
 - Examine wiring for signs of damage or disconnection.
 - **Physical Damage:**
 - Inspect the cabinet for dents or scratches. Contact support if manufacturing defects are identified.
-

7. Product Specifications

- **Material:** Heavy-duty steel for maximum durability.
 - **Dimensions:** 7"D x 14"W x 16"H.
 - **Weight Capacity:** Can hold any portable AED along with a First Aid kit (not included).
 - **Power Requirements:** One 9V battery for lighting system operation.
-

8. Care and Maintenance

- **Routine Cleaning:**
 - Wipe the cabinet with a soft cloth and mild detergent to maintain a clean appearance.
 - Avoid using abrasive materials that could damage the surface.
 - **Regular Inspections:**
 - Periodically check screws, bolts, and seals for signs of wear or loosening.
 - **Security Seals:** Replace the security seals to ensure the cabinet remains tamper-proof.
 - Record inspection details or usage information on the provided inspection cards.
 - **Lighting Maintenance:**
 - Replace the battery in the lighting system annually or as needed to ensure optimal performance.
-

9. Frequently Asked Questions

- **What tools are required for assembly?**
 - A standard screwdriver and drill are recommended.
 - **Can the cabinet be mounted outdoors?**
 - While designed for indoor use, it can be installed in sheltered outdoor areas.
-